

New Gene (Hangzhou) Bioengineering Co., Ltd.

Address: Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District,
Hangzhou, Zhejiang, P.R. China

No. ST2021112801

Statement of the Product Performance on SARS-CoV-2 Mutants

To whom it may concern,

We, **New Gene (Hangzhou) Bioengineering Co., Ltd.** being the official designer and manufacturer of our proprietary **COVID-19 Antigen Detection Kit**, and having factories with address of Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou (310051), Zhejiang, P.R. China, do hereby state that:

The **COVID-19 Antigen Detection Kit** only detects the nucleocapsid antigens of SARS-CoV-2, and the product performance is not affected by mutations in SARS-CoV-2 spike glycoprotein antigens.

The **COVID-19 Antigen Detection Kit** can detect multiple SARS-CoV-2 mutants, including Alpha, Beta, Gamma, Delta, Epsilon, Iota, Kappa, Mu, Fin-796H, and Omicron.

New Gene (Hangzhou) Bioengineering Co., Ltd. reserves all the right to the interpretation of this statement.

New Gene (Hangzhou) Bioengineering Co., Ltd.

Date: 2021/11/28

诺迦(杭州)生物工程有限公司
New Gene (Hangzhou) Bioengineering Co., Ltd.