

**Business Stream Products
Certification Department**



Precisely Right.

TÜV Rheinland LGA Products GmbH · 90431 Nürnberg

Ms. Jinge Li
Zhonghong Pulin Medical
Products Co., Ltd.
West Industrial Park
Luannan County
063500 Tangshan City HEBEI
CHINA

Contact

Tel. +49 911 655-5225
Mail service@de.tuv.com

Date July 26, 2021

Application for : EU-Konformität PSA, MODUL C2
Certificate No. : QP 60157659 Sheet 0001
Device : Protective gloves against chemicals and micro-organisms
according to EN ISO 374-1+A1:2018 and EN ISO 374-5:2016
Type : ZHPFN02
Test requirement : UEReg 425/2016
EN ISO 374-1:2016+A1
EN ISO 374-5:2016

Dear Ms. Li,

The PPE Regulation (EU)2016/425 has fixed rules for the monitoring of manufactured PPE with required checks on the final product described in Annex VII (Module C2) "CONFORMITY TO TYPE BASED ON INTERNAL PRODUCTION CONTROL PLUS SUPERVISED PRODUCT CHECKS AT RANDOM INTERVALS".

The submitted sample of the product has been tested and in this configuration found to be in accordance with the above mentioned requirements.

Enclosed please find the certificate Nr. QP 60157659 0001

This certificate remain valid from 08.03.2021 - 07.03.2024,
date of sample picking: 12.03.2021.

Kind regards

Certification body

Dipl.-Ing. C. Albrecht

Test sample: no, documentation available

TÜV Rheinland
LGA Products GmbH

Tillystraße 2
90431 Nürnberg

Tel. +49 911 655-5225
Fax +49 911 655-5226
Mail service@de.tuv.com
Web www.tuv.com/safety

Board of Management

Dipl.-Ing.
Jörg Mähler, Spokesman

Dipl.-Kfm.
Dr. Jörg Schlösser

Chairman of the
Supervisory Board

Dipl.-Ing.
Ralf Scheller

Nuremberg HRB 26013
VAT No.: DE 811835490

Konformitätsbescheinigung
Gemäß Anhang VII der Verordnung 2016/425/EU
Interne Fertigungskontrolle mit überwachter Produktprüfung

Registrier-Nr.: QP 60157659 0001

Berichts-Nr.: 60221944 006

Hersteller:

Zhonghong Pulin Medical
Products Co., Ltd.
West Industrial Park
Luannan County
Tangshan City
063500 Hebei
P.R. China

Fertigungsstätte:

The Sixth Branch Company of
Zhonghong Pulin Medical Products
Co., Ltd.
Industrial Zone, West Luannan County
Tangshan City
Hebei
P.R. China

Produkt:

Schutzhandschuhe gegen Chemikalien und Mikroorganismen
gemäß EN ISO 374-1+A1:2018 und EN ISO 374-5:2016

Identifikation:

ZHPFN02 Einmalhandschuhe
Farbe: blau
Material: Nitril, Wanddicke ca. 0,05 mm
Größen: S (6,5) - XL (9)
Leistungsparameter: Typ C, Klasse 6, K: NaOH 40%

Dieses Zertifikat dient in Verbindung mit dem o.g. Bericht als Bestätigung, dass eine angemessene statistische Stichprobe der oben genannten hergestellten PSA ausgewählt, beprobt, durch einen stichprobenartigen Labortest gemäß dem o.g. Anhang der EU-Verordnung untersucht, von der Benannten Stelle durchgeführt und bestanden wurde. Somit gewährleistet der Herstellungsprozess die Homogenität der Produktion und die Leistung innerhalb akzeptabler Grenzen.

Gültigkeit:

Von: 2021-03-08 Bis: 2024-03-07

Datum: 2021-07-26

Benannte Stelle



Dipl.-Ing. C. Albrecht

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg

Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-Mail: cert-validity@de.tuv.com <http://www.tuv.com/safety>

Notifiziert unter der Nr. **0197** bei der Kommission der Europäischen Gemeinschaft.

CE Die CE-Kennzeichnung darf bei Einhaltung aller zutreffenden EU-Richtlinien angebracht werden. CE

Certificate of Conformity
According to Annex VII of Regulation 2016/425/EU
Internal production control plus supervised product check

Registration No.: QP 60157659 0001

Report No.: 60221944 006

Manufacturer:

Zhonghong Pulin Medical
Products Co., Ltd.
West Industrial Park
Luannan County
Tangshan City
063500 Hebei
P.R. China

Manufacturing Facility:

The Sixth Branch Company of
Zhonghong Pulin Medical Products
Co., Ltd.
Industrial Zone, West Luannan County
Tangshan City
Hebei
P.R. China

Product:

Protective gloves against chemicals and micro-organisms
according to EN ISO 374-1+A1:2018 and EN ISO 374-5:2016

Identification:

ZHPFN02 disposable gloves
Colour: blue
Material: nitrile, wall thickness ca.0,05 mm
Sizes: S (6,5) - XL (9)
Performance parameter: Type C, class 6, K: NaOH 40%


This certificate in conjunction with above mentioned report serves as confirmation that an adequate statistical sample of the manufactured PPE mentioned above has been selected, sampled, examined by a random basis laboratory test according the above mentioned Annex of the EU Regulation have been carried out by the Notified Body and were passed. Consequently the manufacturing process ensures the homogeneity of production and performs within acceptable limits.

Validity:

from: 2021-03-08 to: 2024-03-07

Date: 2021-07-26

Notified Body



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-Mail: cert-validity@de.tuv.com <http://www.tuv.com/safety>

Notified under No. 0197 to the EC Commission.



The CE marking may be used if all relevant and effective EC Directives are complied with.

